

Nucleic Acid Isolation System  
**QuickGene-810**



# DNA and RNA extraction with a revolutionary 80µm membrane film

QuickGene-810 rapidly extracts DNA/RNA from varied samples with high quality and high yield.

An automated system with extraction kits for reliable results.



Small and self-contained QuickGene-810 takes up minimal space on the lab bench.



Intelligent QuickGene-810 handles samples automatically and extracts DNA or RNA.




FUJIFILM's patented porous membrane makes it possible to extract DNA from eight sets of whole blood samples simultaneously in only six minutes.



FUJIFILM provides several kits for extracting DNA from whole blood and tissue samples and RNA from cultured cells, tissue samples and blood cells.



## FUJIFILM's revolutionary porous membrane (Scanning electron microscope photos)



The system uses a porous, highly adsorptive membrane developed through application of FUJIFILM's advanced polymer membrane production technology. It is only 80µm thick, making it incomparably thinner than conventional glass fibers. QuickGene-810's ultra thin membrane alleviates the risk of contamination from residue in the membrane.


## Easy RNA extraction

Problematic RNA extraction can also be fully automated with QuickGene-810. RNA is much more unstable than DNA, and ribonuclease in the atmosphere or from the operator during the extraction process has sometimes resulted in its degradation. But there is no risk at all of contamination when you use QuickGene-810 because the extraction process occurs automatically in a sealed, enclosed space.

## High purity and high yield without centrifugation


Three pressurizing stages – binding, washing and elution – occur automatically in the unit. Because of the outstanding adsorptive and desorptive properties of the membrane, high-purity nucleic acid can be obtained easily at low pressure without any complex processes such as centrifugation.

### Extraction of DNA from whole blood



### Adsorption of nucleic acids

Owing to their hydrophilic properties, nucleic acids get adsorbed onto the membrane, while proteins and lipids, which are comparatively hydrophobic, tend to seep out of the membrane.



## Processing time (8 samples)

DNA Extraction	
SAMPLE	TIME
WHOLE BLOOD	6 min
TISSUE	13 min


PLASMID DNA Extraction	
SAMPLE	TIME
PLASMID	6 min

RNA Extraction	
SAMPLE	TIME
TISSUE	15 min
CULTURED CELL (adherent / floating)	17 min
CULTURED CELL (6/10cm dish)	11 min
BLOOD CELL	20 min

## High purity

There are almost no impurities in extracted genomic DNA and total RNA. The absence of impurities such as proteins and chaotropic salts means that the extracted products can be used directly in PCR and RT-PCR.

	Purity
DNA	$A_{260/280} > 1.7$
RNA	$A_{260/280} > 1.8$




## High yield

High yields of genomic DNA can be extracted from whole blood and total RNA from cultured cells without any need to use hazardous organic solvents.

### DNA extraction yield compared with competitors

Yield of genomic DNA extracted from 200µL of whole blood (average of ten specimens).



## Simple operation and automated processes reduce risk of contamination

There is no need to set complex extraction conditions. Automation ensures that extraction does not fail because of contamination.



## The specific set of extraction kits supports various samples

QuickGene Extraction kits are optimized for QuickGene system to extract DNA and RNA in the shortest time and with the highest quality. Appropriate kit selectable depending on sample.

Extraction kits			Extraction example
QuickGene DNA whole blood kit S	[ for 96 samples ]	Reference code DB-S	ca.5 µg / Whole blood 200 µl
QuickGene DNA tissue kit S	[ for 96 samples ]	Reference code DT-S	ca.4 µg / 5mg BALB/c Mouse tail
QuickGene Plasmid kit S	[ for 96 samples ]	Reference code PL-S2	ca.12.5 µg / 1ml culture DH5
QuickGene RNA tissue kit S	[ for 96 samples ]	Reference code RT-S2	ca.100 µg / 30mg Mouse liver
QuickGene RNA cultured cell kit S	[ for 96 samples ]	Reference code RC-S	ca.10 µg / 1 × 10 <sup>6</sup> cell HL60 cell
QuickGene RNA cultured cell HC kit S	[ for 96 samples ]	Reference code RC-S2	90 ~ 150 µg / 10cm dish cultured HEK293 cell
QuickGene RNA blood cell kit S	[ for 96 samples ]	Reference code RB-S	ca.4.5 µg / 1 × 10 <sup>7</sup> cell Leukocytes

## Specifications

### Overview

Automated stages : Sample binding, washing and elution  
 Throughput : 1 to 8 samples per run  
 Display : LCD (16 characters × 1 line)

\*Research use only

### Operating conditions

Supply voltage : 100V-240V  
 Power supply frequency : 50/60Hz  
 Operating conditions : Temperature : 15-30  
 Humidity : 30-80% (non-condensing)  
 Power consumption : 65W

### Extraction modes


DNA WHOLE BLOOD  
 DNA TISSUE  
 PLASMID  
 RNA TISSUE  
 RNA TISSUE PLUS  
 RNA CELL  
 RNA CELL PLUS  
 RNA BLOOD

### Extraction kits

QuickGene DNA whole blood kit S ( for 96 samples )  
 QuickGene DNA tissue kit S ( for 96 samples )  
 QuickGene Plasmid kit S ( for 96 samples )  
 QuickGene RNA tissue kit S ( for 96 samples )  
 QuickGene RNA cultured cell kit S ( for 96 samples )  
 QuickGene RNA cultured cell HC kit S ( for 96 samples )  
 QuickGene RNA blood cell kit S ( for 96 samples )

### Physical specifications

Dimensions : 450 (w) × 330 (D) × 400 (h) mm  
 Weight : 21kg



# FUJIFILM

FUJIFILM Corporation 7-3, Akasaka 9-Chome Minato-ku, Tokyo 107-0052, Japan, Tel: +81-3-6271-2158, Fax: +81-3-6271-3136, E-mail: sginfo@fujifilm.co.jp

FUJIFILM Europe GmbH Heesenstr. 31, 40549 Dusseldorf, Germany, Tel: +49-211-5089-174, Fax: +49-211-5089-9144, E-mail: lifescience@fujifilm-europe.de

FUJIFILM UK Ltd. Unit 12 St Martins Way, St Martins Business Centre, Bedford, MK42 0LF, U.K., Tel: +44-1234-572226, Fax: +44-1234-245293, E-mail: lifesciences@fuji.co.uk

富士胶片(中国)投资有限公司 31st floor, Hong Kong New World Tower, No.300 Huai Hai Zhong Road, Shanghai, P.R China, Tel: +86-21-3302-4655 ext.363, Fax: +86-21-6384-3322, E-mail: wgxiang@fujifilm.com.cn

<http://lifescience.fujifilm.com>

Specifications and system configuration subject to change for improvement without notice. All other product names mentioned herein are the trademarks of their respective owners.

